

# Spectrometry\_Gamma

## A New Module For The EduGATE Project

Uwe Pietrzyk & Michaela Gaens

Institute of Neurosciences and Medicine (INM), Research Center Juelich, Germany  
Department of Mathematics and Natural Sciences, University of Wuppertal, Germany

Contact: [u.pietrzyk@fz-juelich.de](mailto:u.pietrzyk@fz-juelich.de) or [uwe.pietrzyk@uni-wuppertal.de](mailto:uwe.pietrzyk@uni-wuppertal.de)

(Basic Paper → „Zeitschrift für Medizinische Physik“ / Z. Med. Phys. 23 (2013) 65-70)


- *Long is the way of theory,  
short and effective by examples –*

(Lucius Annaeus Seneca (the Younger),  
Epistulae Morales)

- **Spectrometry\_Gamma:**


- introduces to the basics of gamma spectra obtained from a scintillation detector
- Inspired by **chapter 10** of  
***Physics in Nuclear Medicine***  
**(Cherry, Sorensen & Phelps, 2012)**
- provides GATE-macros
- comes in two versions: **Basic & Extended**
- Programme Code in C and for ROOT to calculate and plot relevant parameters during **Compton Scatter**

# Pulse-height spectrum for $^{137}\text{Cs}$ showing typical structures



**FIGURE 10-3** Actual pulse-height spectrum recorded with a NaI(Tl) detector and  $^{137}\text{Cs}$  (662-keV  $\gamma$  rays, ~30 keV Ba x rays). Compare with Figure 10-2B.

from: Cherry et al, 2012



(Cherry, et al, 2012)

## Basics of Compton\_Scatter

$$E_{sc} = E_0 / [1 + (E_0 / 0.511) (1 - \cos\theta)]$$

$$E_{re} = E_0 - E_{sc}$$

for 180-degree scattering:

$$E_{sc}^{\min} = E_0 / [1 + (2 E_0 / 0.511)]$$

$$E_{re}^{\max} = E_0 - E_{sc}^{\min}$$

# running: root -l compton\_scatter\_root.c

```
for (int i=0; i<n_values; i++)
{
// energy of the incident gamma
E_gamma[i] = (i+1) * e_step;

// energy of the recoil electron; -> "compton electron" -> "compton edge"
Emax_re[i] = E_gamma[i] * E_gamma[i] / (E_gamma[i] + m_electron_2);

// energy of the back-scattered gamma
Emin_sc[i] = E_gamma[i] - Emax_re[i];

printf( " E_gamma: %f Emax_re: %f Emin_sc: %f (MeV)\n", E_gamma[i], Emax_re[i], Emin_sc[i]);

}
```

--> output


```
E_gamma: 0.010000 Emax_re: 0.000377 Emin_sc: 0.009623 (MeV)
E_gamma: 0.020000 Emax_re: 0.001452 Emin_sc: 0.018548 (MeV)
E_gamma: 0.030000 Emax_re: 0.003152 Emin_sc: 0.026848 (MeV)
E_gamma: 0.040000 Emax_re: 0.005415 Emin_sc: 0.034585 (MeV)

.
.


E_gamma: 0.660000 Emax_re: 0.475806 Emin_sc: 0.184194 (MeV)
E_gamma: 0.670000 Emax_re: 0.485035 Emin_sc: 0.184965 (MeV)
E_gamma: 0.680000 Emax_re: 0.494281 Emin_sc: 0.185719 (MeV)
E_gamma: 0.690000 Emax_re: 0.503543 Emin_sc: 0.186457 (MeV)
E_gamma: 0.700000 Emax_re: 0.512820 Emin_sc: 0.187180 (MeV)
```

**Plot created via: root -l compton\_scatter\_root.c**


## Compton scattering at 180-degree



## Compton scattering at 180-degree




# Pulse-height spectrum for $^{137}\text{Cs}$ showing typical structures




**FIGURE 10-3** Actual pulse-height spectrum recorded with a NaI(Tl) detector and  $^{137}\text{Cs}$  (662-keV  $\gamma$  rays, ~30 keV Ba x rays). Compare with Figure 10-2B.

from: Cherry et al, 2012

# Spectrometry\_Gamma\_Basic




# Basic Setup



- Two opposing NaI-crystals allow to detect:
  - back-scattered gamma
  - characteristic X-ray, escaping from the crystal (head 0)


*simplified* Gamma Camera → 2 **crystals**, no collimator,  
**collimated gamma source** (rays directed towards head 0 !!)

# Compton scattering at 180-degree




E\_gamma: 0.510000 Emax\_re: 0.339778 Emin\_sc: 0.170222 (MeV)

**511 keV**



**note: different scales!!**




# Pulse-height spectrum for

$^{197}\text{Hg}$




some more details!!



**FIGURE 10-4** Pulse-height spectrum for  $^{197}\text{Hg}$  ( $E_{\gamma} = 77.3 \text{ keV}$ ) recorded with NaI(Tl). Iodine escape peak (45-50 keV) is due to escape of characteristic iodine x ray (~30 keV) following a photoelectric absorption event in detector.


**80 keV**




**note: different scales!!**

**escape peak:  $E_0 - E_{\text{X-rays}}$**

**Iodine (in NaI-crystal):**  
**characteristic X-rays**  
**following**  
**photoelectric**  
**absorption**




**$E_{\text{gamma}}: 0.080 \quad E_{\text{max\_re}}: 0.019076 \quad E_{\text{min\_sc}}: 0.060924 \quad (\text{MeV})$**



[http://xdb.lbl.gov/Section1/Periodic\\_Table/I\\_Web\\_data.htm](http://xdb.lbl.gov/Section1/Periodic_Table/I_Web_data.htm)

K and L shell emission lines in electron volts

Ka <sub>1</sub>	Ka <sub>2</sub>	Kb <sub>1</sub>	La <sub>1</sub>	La <sub>2</sub>	Lb <sub>1</sub>	Lb <sub>2</sub>	Lg <sub>1</sub>
28,612.0	28,317.2	32,294.7	3,937.65	3,926.04	4,220.72	4,507.5	4,800.9



E\_gamma: 0.080 Emax\_re: 0.019076 Emin\_sc: 0.060924 (MeV)

# Starting Basic Setup (Gate\_v7.0)

in a shell-window type:  
config\_starter\_mac\_70.sh

EduGate Spectrometry\_Gamma Configuration

CameraType	2heads
ViewPointThetaPhi	0 90
VisuOnOff	novisu
SourceActivity	100 Bq
CrystalMaterial	Nal
SourceVolMaterial	Vacuum
SourceType	gamma_collim
SourceEnergy	80
E_res	0.10

[Generate configuration.mac and Start](#)

specify:

energy &

energy resolution

# Starting Basic Setup (Gate\_v7.0)

can change selection presets in **Spectrometry\_Gamma.txt**



```
CameraType: 2heads; 1head; 4heads;  
ViewPointThetaPhi: 0 90; 90 0; -90 0; 89 90; 15 30; 30 30; 45 45; 60 60;  
VisuOnOff: novisu; visu;  
SourceActivity: 100 Bq; 1000. Bq; 10000. Bq; 100000. Bq;  
CrystalMaterial: NaI; BGO; LSO; GSO; PWO; LuAP; YAP; CZT; Silicon; Lead; Tungsten;  
SourceVolMaterial: Vacuum; Water; Plexiglass; Air; PVC;  
SourceType: gamma_collim; gamma;  
SourceEnergy: 80; 10; 50; 60; 90; 100; 120; 140; 160; 200; 240; 364; 511; 1000; 1600;  
E_res: 0.10; 0.00; 0.01; 0.05; 0.15; 0.20; 0.25;
```

# Result: Basic Setup (Gate\_v7.0)

