




EduGATE → follows the basic Idea of a Development / Processing Pipeline ...





Signal Processing Pipeline: Sensitive Detectors / Digitizer



- Pre-programmed components
 - Sensitive detectors
 - Trajectory analyser
- Digitizer
 - Linear signal processing chain
 - Modular: set-up via scripting


```
/gate/digitizer/modules/insert blurring  
/gate/digitizer/blurring/setResolution 0.15  
  
/gate/digitizer/modules/insert spblurring  
/gate/digitizer/spblurring/setResolution 0.15  
  
/gate/digitizer/modules/insert readout  
  
/gate/digitizer/modules/insert thresholder  
/gate/digitizer/thresholder/setThreshold 50. keV
```

Basic GATE Macro → Phantom and Source

PHANTOM VOLUME

Create the attenuating volume/**phantom**

```
/gate/world/daughters/name Phantom  
/gate/world/daughters/insert cylinder  
/gate/Phantom/geometry/setRmax 50.0 mm  
/gate/Phantom/geometry/setRmin 49.0 mm  
/gate/Phantom/geometry/setHeight 20. cm  
/gate/Phantom/placement/setTranslation  
          0. 0. 0. cm  
/gate/Phantom/setMaterial Water  
/gate/Phantom/vis/setColor blue  
#/gate/Phantom/vis/forceWireframe  
/gate/Phantom/vis/forceSolid  
  
/gate/Phantom/attachPhantomSD
```



SOURCE Volume

```
# Add an extra object for source volume  
/  
gate/world/daughters/name source_vol  
/gate/world/daughters/insert cylinder  
/gate/source_vol/geometry/setRmax 2. cm  
/gate/source_vol/geometry/setRmin 0. cm  
/gate/source_vol/geometry/setHeight 8. cm  
/gate/source_vol/placement/setTranslation  
          0. 0. 0. cm  
/gate/source_vol/setMaterial Water  
/gate/source_vol/vis/setColor magenta  
/gate/source_vol/vis/forceWireframe  
#/gate/source_vol/vis/forceSolid  
  
/gate/source_vol/attachPhantomSD
```

file: Gamma_Camera.txt
with start-settings

→ large number of combinations

ViewPointThetaPhi: 0 90; 90 0; -90 0; 89 90;
 VisuOnOff: novisu; visu;
 Src_Act: 20 Bq; 50 Bq; 100 Bq; 1000 Bq; 10000 Bq; 0.1 MBq; 0.5
 CameraType: camera_Tc; camera_I_131;
 x_pos: 20.0; 15.0; 10.0;
 CollMat: Lead; Vacuum; Air; Copper; Iron; Tungsten; Plexiglass;
 CrysMat: NaI; BGO; LSO; GSO; PWO; LuAP; YAP; CZT;
 PhanMat: Plexiglass; Water; Air; Vacuum; Lead; PVC; Copper;....;
 PhanRmax: 50 mm; 52 ; 54 ;
 PhanRmin: 49 mm; 49.99 mm; 49.9 mm; 49.5 mm;
 SrcVolMat: Plexiglass; Air; Water; Vacuum; PVC;
 SrcType: src_gamma_3_lim_ang; src_gamma_lim_ang;
 Src_E: 140; 80; 364; 511; 1000; 2284;
 E.blur: 0.1; 0.0; 0.05; 0; 0.15; 0.20; 0.25;
 SP.blur: 2.0; 0.5; 1.0; 1.5; 2.5; 3.0; 4.0; 5.0;
 E_low: 5; 50; 100;
 E_up: 2000; 100;
 PileupOnOff: no_pile_up; pile_up;
 t_slice: 10.; 1.0; 0.1; 0.01; 0.001; 0.0001;
 t_stop: 10.; 1.0; 0.1; 0.01; 0.001; 0.0001;

default selections


EduGATE Configuration Menu

ViewPointThetaPhi	0 90
VisuOnOff	novisu
Src_Act	20 Bq
CameraType	camera_Tc
x_pos	20.0
CollMat	Lead
CrysMat	NaI
PhanMat	Plexiglass
PhanRmax	50 mm
PhanRmin	49 mm
SrcVolMat	Plexiglass
SrcType	src_gamma_3_lim_ang
Src_E	140
E.blur	0.1
SP.blur	2.0
E_low	5
E_up	2000
PileupOnOff	no_pile_up
t_slice	10.
t_stop	10.

file: configuration.mac

← generated based on selections

```
/control/alias ViewPointThetaPhi 0 90
/control/alias VisuOnOff novisu
/control/alias Src_Act 10000 Bq
/control/alias CameraType camera_Tc
/control/alias x_pos 20.0
/control/alias CollMat Lead
/control/alias CrysMat NaI
/control/alias PhanMat Plexiglass
/control/alias PhanRmax 50 mm
/control/alias PhanRmin 49 mm
/control/alias SrcVolMat Plexiglass
/control/alias SrcType src_electron_lim_ang
/control/alias Src_E 140
/control/alias E.blur 0.1
/control/alias SP.blur 2.0
/control/alias E_low 5
/control/alias E_up 2000
/control/alias PileupOnOff no_pile_up
/control/alias t_slice 10.
/control/alias t_stop 10.
```



Hints for bookkeeping

construct Root-file-name from specific parameters:

`GC_camera_Tc_Lead_Nal_Plexiglass_Plexiglass_src_gamma_3_lim_ang_140_1000Bq_no_pile_up_0.1_2.0_5_2000_1.0_10..root`

/control/alias RootFileName GC_{CamType}_{CollMat}_{CrysMat}_{PhanMat}_{SrcVolMat}_{SrcType}_{Src_E}_{Src_Act}
{PileupOnOff}{E.blur}_{SP.blur}_{E.low}_{E.up}_{t.slice}_{t.stop}


GC → Gamma Camera:

(Correspondence - selected items)

CamType	camera_Tc
SrcType	sources_gamma_3
CollMat	Lead
CrysMat	Nal
Src_E	140 keV
Src_A	0.1MBq
E.blur	0.0
SP.blur	2.0 mm
E_low	5 keV
E_up	2000 keV
t_slice	1.0 sec → # runs
t_sop	10 sec → simulation time

Basic EduGATE Examples


Planar Gamma Camera System



- Left: Magenta: cylindrical Phantom & Source, Blue: hollow cylinder serving as attenuating material of various types,
- Right: Head – Blue: Collimator – Yellow: Crystal, Grey: Back-compartment, Table
 - Green: Gamma-rays


Basic EduGATE Examples

Planar Gamma Camera System



Perspective View – Facing the Front: Magenta: cylindrical Phantom & Source,
Blue: hollow cylinder serving as attenuating material of various types,
White: Camera-Head, Green: Gamma-rays

Total Spectrum of the detected events



GC_src_gamma_3_lim_ang_Lead_NaI_140_100 MBq_pile_up_0.0_2.0_5_2000_0.1_1.0.root

additional line for the DIGITIZER:

/gate/digitizer/Singles/insert pileup

/gate/digitizer/Singles/pileup/setDepth 1

/gate/digitizer/Singles/pileup/setPileup 100 ns

no energy blurring

Demo-1(Gamma Camera)

- How to start??
- → use basic scripts to setup environment:

```
source config_starter_80_mac.csh
```

```
#!/bin/csh  
# file config_starter80.csh (also available for GATE 7.2)
```

```
setenv CURPATH $cwd
```

```
cd /Applications/root_v6.08.06/bin  
source thisroot.csh
```

```
cd /Applications/geant4.10.03.p01-install/bin  
source geant4.csh
```

```
setenv PATH /Applications/gate_v8.0-install/bin:$PATH
```

```
cd $CURPATH
```

```
# select configuration:
```

```
root -l 'GenerateGateConfiguration.C( "Gamma_Camera.txt" )'
```


```
# start GATE:
```

```
Gate Gamma_Camera.mac | & tee Gamma_Camera.log
```

```
echo Time is: `date` | & tee -a Gamma_Camera.log
```

```
# start Analysis:
```

```
root -l Gamma_Camera.C → Select RootFile ...
```



```
##### Number of emitted particles : 300438
##### Number of detected events : 1428
##### Primary events : 1185
##### Scatter in the phantom : 42
##### Scatter in the table : 0
##### Scatter in the collimator : 7
##### Scatter in the crystal : 168
##### Scatter in the backcompartment : 25
```

Contact

Prof. Dr. Uwe Pietrzyk

Institute of Neuroscience and Medicine (INM-4)

- Medical Imaging Physics -

Forschungszentrum Jülich GmbH

U.Pietrzyk@fz-juelich.de

<http://www.fz-juelich.de>

and

School of Mathematics and Natural Science (C)

- Medical Physics -

Bergische Universität Wuppertal

Uwe.Pietrzyk@uni-wuppertal.de

<http://www.medizinphysik.uni-wuppertal.de>

Germany